



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
TRAINING STAIRS

Product Code :
SAMCOMEDICAL210003A015



Description :

TRAINING STAIRS Sturdily built in two sections for installation in a straight single side. Steps arrangement : Four 6" (15cm.) steps lead to a platform size 30" x 30" x 24" high. Steps. : 30" wide and 11" deep. Hand Rails. : All-around hardwood handrails at two different Height to accommodate adults and children. Construction. : Steps and platform made of commercial board, hand rails with supporting bars made out of hardwood. Polished natural wood finish



Samco Medical

Website: www.samcomedical.com,