



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
CURBS AND RAMP

**Product Code :**  
SAMCOMEDICAL210003A016



**Description :**

CURBS AND RAMP Hardwood constructed training set consists of ONE ramp and THREE nested curbs stored below the ramp. The ramp graduates from 20mm to 200mm and is 84" X 44". The curbs are 40"x 40"x 6", 36"x 36"x 4" and 32"x 32"x 2" respectively. Natural wood finish with anti-slip safety mats on ramp and curbs without handrails



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),