



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
STATIC CYCLE EXERCISER

Product Code :
SAMCOMEDICAL210003A017



Description :

STATIC CYCLE EXERCISER Ruggedly built tubular steel frame stands solidly on four rubber tipped properly balanced legs. Provided with adjustable standard cycle seat and specially designed adjustable handle. Fitted with one hard rubber tyred cycle wheel and standard chain and sprocket for bicycle drive. Comes with chrome plated wing nut, pressure rod and ball bearing resistance roller which permits controlled variation in riding effort. Odometer shows speed and kilometer done. Oven baked finish with working parts finish with chrome plating



Samco Medical

Website: www.samcomedical.com,