



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
CYCLE ERGOCISER, Semi-Recumbent

Product Code :
SAMCOMEDICAL210003A019



Description :

CYCLE ERGOCISER, Semi-Recumbent The adjustable semi-recumbent ergociser provides a biomechanical seat for those patients requiring back support & reduced gluteal pressure. The sturdy, heavy duty frame with powder coated finish and non-friction magnetic resistance system provides smoothness, durability and runs virtually silent. FEATURES: Frame: Heavy gauge pipe frame with heavy fly wheel has Powder coated finish. Resistance: Non-friction, Magnetic, with EIGHT levels of resistance system. Display: Time, Distance, Speed & Calories



Samco Medical

Website: www.samcomedical.com,