



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
ARM CONTRAST BATH

Product Code :
SAMCOMEDICAL210004A001

Description :

ARM CONTRAST BATH The double-walled, well insulated, Stainless steel sheet constructed bath is highly suitable for treatment of the Arms. The bath unit has Two tubs mounted on a wide base and fitted with waterinlet and over-flow pipes. Works on 220V AC. Overall Size : 90cm x 80cm x 95cm



Samco Medical

Website: www.samcomedical.com,