



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
COLD THERAPY UNIT

Product Code :
SAMCOMEDICAL210004A002



Description :

COLD THERAPY UNIT Made out of heavy gauge stainless steel sheet. Tank is fully insulated in between top and sides. Fitted with heavy duty cooling unit and thermostat for temperature control. A drain valve is provided for emptying. Mounted on 3" sizes ball bearing rubber castors for easy mobility. Comes with 6 standard packs (Size 25 x 30 cm), containing special silica gel crystals. Finished bright inside and special paint outside



Samco Medical

Website: www.samcomedical.com,