



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
CONTRAST BATH

Product Code :
SAMCOMEDICAL210004A003



Description :

CONTRAST BATH Suitable for treatment of Arms & Legs. Two Stainless Steel heavy gauge tanks with inside measurements of 50cm x 35cm x 38cm deep each are mounted on 10 cm size ball bearing rubber castors. One tank is fitted with special heating unit and the other with heavy duty cooling unit. Each tank provided with Digital Thermometer and Thermostat for Temperature Control, Pilot Lights & Water Outlet. Both tanks are well insulated and placed side by side in a heavy tubular steel frame covered with painted steel sheet. Works on 220 V. A.C. UNIT COMES WITH 5 kv. VOLTAGE STABILISER



Samco Medical

Website: www.samcomedical.com,