



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
INFRA RED LAMP

Product Code :
SAMCOMEDICAL210004A004



Description :

INFRA RED LAMP Floor stand model for IR radiation treatments. IR bulb of 150w is fitted in metal casing with auto-matic cut off timer of 15 minutes. Adjustable arm permits wide range of adjustments. Mounted on wide base with ball bearing casters. Comes with suitable power cord



Samco Medical

Website: www.samcomedical.com,