



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
ULTRA VIOLET LAMP

Product Code :
SAMCOMEDICAL210004A005



Description :

ULTRA VIOLET LAMP Floor stand model for Ultraviolet radiation treatments. UV bulb of 300w is fitted in metal casing with auto-matic cut off timer of 15 minutes. Adjustable arm permits wide range of adjustments. Mounted on heavy duty base fitted with ball bearing casters. Comes with suitable power cord



Samco Medical

Website: www.samcomedical.com,