



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
MOIST HEAT THERAPY UNIT

Product Code :
SAMCOMEDICAL210004A006



Description :

MOIST HEAT THERAPY UNIT Made out of heavy gauge stainless steel sheet double walled well insulated in between. Overall size 38cmx 30cmx 45cm high. Fitted with 1000 watts immersion heater, pilot light and thermostat for heat control. Comes with 4 standard size 25cmx30cm steam packs containing special heat retaining silica gel crystals. Polished finish inside and special paint finish outside



Samco Medical

Website: www.samcomedical.com,