



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
SUSPENSION FRAME

**Product Code :**  
SAMCOMEDICAL210006A001

**Description :**

SUSPENSION FRAME Tubular steel construction, overall sizes 210cm long x 110cm wide at the base, 50cm wide at the top and 190cm high. All welded steel rod mesh at the top for complete suspension with slings & ropes, 12 hooks are provided on each side to hold suspension gear. Mesh facilitate suspension at any desired point in relation to the patient and as required in each individual case. Comes complete with Suspension Gear comprising 9 canvas slings (6 limb, 2 trunk & 1 head), 8 single ropes & four double pulley ropes with necessary hooks & cleats, 18 graduated springs of different tension, 12 S'hooks and 2 ankle straps. Spray enamel finish



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),