



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
CERVICAL & LUMBER TRACTION PACKAGE

Product Code :
SAMCOMEDICAL210006A003



Description :

CERVICAL & LUMBER TRACTION PACKAGE Microprocessor controlled Traction system, (8) Eight Traction Modes and operation through a large LCD display & Programmable Memory but with Indigenous NML-2731 ELECTRIC HI_LO RACTION TABLE. TRACTION MACHINE FEATURES : Traction Force: 90 kgs. Maximum. Traction Modes (8): Intermittent, Static, Harmonic Progressive, Regressive. Traction Speed: Variable Traction Speed. Traction Memory: Stores 20 Treatment Patterns. Traction Safety: A Self- diagnosis functions. : Checking of the system when high force is set above 10 kg. : Automatic tension release when emergencies occur. : Display of error contents and solutions. SPECIFICATIONS : Traction Force: Hold 1.8-90 kg. Hold- Rest Time: 1-99 sec. Independently set. Treatment Time: 1 - 99 minutes. Traction Speed: 9 - 35 seconds to reach 90 kg. Zero Calibration: Automatic Safety Class: Class 1 Type BF, IEC 601. TRACTION TABLE with VARIABLE HEIGHT (Indigenous) : Electrically Operated, Three section Hi-Low Table. Head End.: With Nose Cut-out, and Adjusts from -20 to + 20degree. Lumber Section.: Friction Free Rolling Lumber Section with adjustment Machine Mount Top. Height: Electrically Adjusts from 55cm to 100cm. Accessories.: Pelvic & Thoracic belt set with Bar



Samco Medical

Website: www.samcomedical.com,