



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
AKAY HI-TRAC UNIT

**Product Code :**  
SAMCOMEDICAL210006A004



**Description :**

AKAY HI-TRAC UNIT Akay Hi-Trac micro-computer controlled unit for INTERMITTENT/CONSTANT, LUMBER-CUM-CERVICAL TRACTION is based on solid-state circuitry. Three section, Fixed Height Bed with storage cabin underneath is an integral part of the unit. TECHNICAL SPECIFICATION: Therapy Modes: 8 Traction Modes for either Lumber or Cervical Traction, in Progressive, Regressive Intermittent and Static etc. Memory/Programme: 30 Treatment Patterns can be programmed & Stored in non volatile memory. (24 Programs are stored and Six programs are free). Traction Force: Adjustable from 5 kg to 60 kg. Base Force: Adjustable from 0 kg to 54 kg. Traction Hold Time: Adjustable from 0 to 99 Sec. Base Trac Hold Time: Adjustable from 0 to 99 Sec. Treatment Time: Continuously adjustable from 0 to 99 minutes, with acoustic signal & automatic reduction of traction force. Displays: Digital display of Traction Force, Base Force, TF Hold Time BF Hold Time, Treatment Time & Graphical Display of Trac. Mode. Accessories: Adjustable Flexion Stool, Head-Halter, Pelvic & Thoracic Belts. Couch/Bed: Fixed Height-70mm, foam padded top is divided in Three sections. Head section has an adjustable incline. Auto-Tension: The set traction force is Adjustment automatically maintained throughout operation. Patient Safety: Patient Stop, hand held switch. Power Supply: 220-240 V AC.(use stabilizer



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),