



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
NICE-TRAC TRACTION UNIT

Product Code :
SAMCOMEDICAL210006A005



Description :

NICE-TRAC TRACTION UNIT An Intermittent traction machine is an accurate & sturdy unit. Works on solid state circuitry with reliable, silent & smooth geared motor. SPECIFICATIONS: Traction Force: Direct Pulling 5 to 45 kgs. With Doubler upto 90 kgs. Hold Time: 20, 40, 60 & 80 seconds. Rest Time: 1, 5, 10 & 20 seconds. Treatment Time: 60 mins. with digital display. Alarms: At the end of treatment time. Safety: Patient safety switch. Mounting: Table or Wall mount. Accessories: Cervical Head Holder with Bar, Lumber Traction Belts with Bar, Main Cord & Pulley Doubler



Samco Medical

Website: www.samcomedical.com,