



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
TRACTION SET CERVICAL

Product Code :
SAMCOMEDICAL210006A007



Description :

TRACTION SET CERVICAL Comfortable Hi-Lo Electric, THREE section table is suitable for TRACTION & TREATMENT Manipulation. Suitable for NML-2721 OR NML-2725 Traction Unit. FEATURES: Height Adjustment: Electrically, with hand switch, from 50cm to 95cm. Upholstered Top Size: 195cm long x 70cm wide. Sections: THREE Sections, drop end section with Nose cut out adjusts from -20 to + 20 deg. Non-friction rolling Lumber section Tilts 60 degree above horizontal. Traction Support: Traction unit can be bolted to an end, for Cervical & Lumber traction. Accessories: Traction Belts, Head-Halter and Flexion stool. Power Supply: 220V AC



Samco Medical

Website: www.samcomedical.com,