



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
HI-LOW MAT PLATFORM

Product Code :
SAMCOMEDICAL210007A004



Description :

HI-LOW MAT PLATFORM This HI-LOW Mat-Platform is ideal for range of motion, transfer, strength building and other active therapies. Platform with height adjustment and adjustable backrest makes delivering therapy easier on the therapist. SPECIFICATIONS: Frame-Work: Square tube constructed, with wide base and finished oven baked enamel. Height Adjustment: Electrically Adjustable with a Hand switch from 51cm to 95cm. Mat Platform: Foam padded, with adjustable Back-rest, Size 120cm x 210cm. Motor System: Imported with vibration free movement. Weight capacity 120 kgs



Samco Medical

Website: www.samcomedical.com,