



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
TILT TABLE

Product Code :
SAMCOMEDICAL210007A005



Description :

TILT TABLE Electrically operated Tilt Table with foam padded top, provided with three straps to hold the patient. Heavy duty motor provides silent and smooth lifting action with push button control switches. Range of tilt is calibrated from 0 - 90 degrees. The table top is 61cm wide x 195cm long x 80cm high, and is fitted on a heavy duty tubular steel frame with four locking wheels of 10 cm dia. for easy mobility. Table comes complete with an Utility Tray and gripping handle for ADL activity. Full 90 degree of tilt plus negative facility also Steel frame oven baked finish. Works on 220 V AC



Samco Medical

Website: www.samcomedical.com,