



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Continuous Passive Motion Unit

Product Code :
SAMCOMEDICAL210007A008



Description :

Continuous Passive Motion Unit

Technical Specification :

CONTINUOUS PASSIVE MOTION UNIT The CPM is designed for upper extremities (Elbow), to stimulate the healing of cartilage in addition to the reduction of post operative swelling & pain. The solid state circuitry control panel with highly visual LED digital readouts, allows the Flexing & Extending an Elbow joint through a pre-fixed range of Motion & Time. Wrist rotation and Wrist roll movement can also be achieved. Unit is fitted on mobile Electrical Jack to adjust height according to patient's requirement. The necessary limb supports provide a comfortable surface for prolonged contact. Comes with necessary accessories for Supination & Pronation, Wrist Roll and Elbow passive motions



Samco Medical

Website: www.samcomedical.com,