



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
VARUN CONTINUOUS PASSIVE MOTION UNIT

Product Code :
SAMCOMEDICAL210007A009



Description :

VARUN CONTINUOUS PASSIVE MOTION UNIT The CPM is designed for lower extremities (Knee & Hip), to stimulate the healing of cartilage in addition to the reduction of post operative swelling and pain. The solid state circuitry control panel with highly visual LED Digital readouts allows Flexing and Extending a joint through a pre-fixed range of motion and time. Proper fit of Calf and Thigh lengths are easily achieved. The necessary limb support provides a comfortable surface for prolonged contact and are washable in cold water. The light weight unit is constructed of Aluminium Sheet and is fitted with carrying Handle and Bed fixing clamps.

SPECIFICATIONS: Flexion Angle: Digital display of Flexion & Extension Angles. 0 to 108 deg. Timer: Digital Timer upto 99 minutes. Angle Adjustment: Two Inching knobs to adjust starting angle. Speed: Fixed single speed. Right or Left Leg can be exercised by little adjustments. Patient Safety: Patient safety switch provided



Samco Medical

Website: www.samcomedical.com,