

**Product Name :**  
Samcomedical Continuous Passive Motion Unit

**Product Code :**  
SAMCOMEDICAL210007A010



**Description :**

Samcomedical Continuous Passive Motion Unit

**Technical Specification :**

Clone of CONTINUOUS PASSIVE MOTION UNIT The CPM is designed for upper extremities (Elbow ), to stimulate the healing of cartilage in addition to the reduction of post operative swelling & pain. The solid state circuitry control panel with highly visual LED digital readouts, allows the Flexing & Extending an Elbow joint through a pre-fixed range of Motion & Time. Wrist rotation and Wrist roll movement can also be achieved. Unit is fitted on mobile Electrical Jack to adjust height according to patient's requirement. The necessary limb supports provide a comfortable surface for prolonged contact. Comes with necessary accessories for Supination & Pronation, Wrist Roll and Elbow passive motions.

**SPECIFICATIONS:**

Flexion Angles: Adjustable from 0 - 160 degree & displayed by Bargraph. Inching: Flexion or Extension Starting position can be set by switches. Timer: Digitally displayed, 0 - 99 minutes. Height Adjustment: Electrically, by Push button switches. Right or Left arm can be exercised by little adjustments. Safety: Patient hand held safety switch

---

Website: [www.samcomedical.com](http://www.samcomedical.com),

