



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
MULTI EXERCISE THERAPY UNIT

Product Code :
SAMCOMEDICAL210003A028



Description :

MULTI EXERCISE THERAPY UNIT Bright polished Teakwood construction,suitable for Shoulder, Elbow, Chest and lower limb exercises. Unit Consists of the following attachments. SPECIFICATIONS: Wall bars 240 cm high x 75 cm wide. Detachable wrist roller with side handle for supination and pronation exercises. Removeable rowing machine cum sliding seat. Unit is fitted with 12 pulleys and Two sets of locking weights (10, one kilo plates in each set) that allows effective weight resistive exercises at different levels



Samco Medical

Website: www.samcomedical.com,