



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
EXERCISE CHAIR

**Product Code :**  
SAMCOMEDICAL210003A029



**Description :**

EXERCISE CHAIR Specifically designed chair, which require little space and can be used in small therapy departments where space is limited. Tubular steel constructed chair of 100 cm long x 67 cm wide is fitted with overhead bar that extends out 80cm from the frame. The chair with two stabilizing straps, three handles two foot straps, ropes with pulleys and overhead pulley system can be used for variety of exercises involving Reciprocal Motion for Elbow, Shoulder, Knees and Hips; Active Quadriceps and Hamstrings Stretching; Resistive and Assistive exercises. The set of Ten locking weights provides graded resistance for all exercises. Spray enamel finish



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),