



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
AXIAL SHOULDER EXERCISER

Product Code :
SAMCOMEDICAL210003A031



Description :

AXIAL SHOULDER EXERCISER For shoulder & supination-pronation exercises. Used in the same way as a conventional shoulder wheel & wrist machine. A wall mounting unit suitable for small clinics is adjustable in height & radius and fitted with calibrated sensitive resistance mechanism, resistance is controllable from zero to maximum by tightening the strap. A 360 deg. scale is provided on the drum to measure and record the degree of revolution from either direction. Adjustment bars are chrome plated. Comes with wall fixing plates



Samco Medical

Website: www.samcomedical.com,