



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
ROWING MACHINE cum SLIDING SEAT

Product Code :
SAMCOMEDICAL210003A036



Description :

ROWING MACHINE cum SLIDING SEAT Tubular steel frame, bright chrome plated finish, size 60"long x 13"wide. Specifically designed foam padded seat is fitted with FOUR nylon castors and slides smoothly over the frame. Comes with FIVE Springs of different tensions. TWO high tension springs attached to the seat for leg exercises and Three light tension springs are fitted with Two Handles for Rowing exercises. Resistance can be varied by removal or adding the springs



Samco Medical

Website: www.samcomedical.com,