



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
ROTARY WRIST MACHINE

**Product Code :**  
SAMCOMEDICAL210003A039



**Description :**

ROTARY WRIST MACHINE Single grip unit of precision machined metal drum with calibrated sensitive resistance control mechanism, is fitted on a height adjustable bar to suit different patients. Resistance can be increased or decreased for Supination and Pronation exercises. Range of motion can also be measured and recorded from the scale provided on the drum



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),