



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
WRIST ROLL

**Product Code :**  
SAMCOMEDICAL210003A040



**Description :**

WRIST ROLL One piece aluminum roll in three different diameters i.e. 30mm, 40mm and 50mm fitted to wall board size: 20cm x 90cm long with adjustable friction control mechanism, controllable from zero to maximum for flexion and extension exercises of wrist and forearm. Oven baked finish. Wooden parts finished natural



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),