



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
WRIST CIRCUMDUCTOR

**Product Code :**  
SAMCOMEDICAL210003A041



**Description :**

WRIST CIRCUMDUCTOR Cast iron wheel of 30cm diameter with drum stand is mounted on a laminated wooden base size:65cm x 22cm. Has a small adjustable handle fitted to the wheel. Knob is provided on the drum to control resistance from zero to maximum. A hollowed foam padded platform size: 22cm x 13cm x 15cm high is fitted with two straps to hold arm and is fitted over the base in front of the wheel for wrist circumduction. Metal parts oven baked finish



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),