



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
GRIP EXERCISER

Product Code :
SAMCOMEDICAL210003A042



Description :

GRIP EXERCISER Chrome plated finished, metallic grip exerciser is fitted on a 50mm wide laminated base. Suitable for patients with very weak / infirm hand. Gives heavy resistance with Six detachable springs. The upper bar is calibrated in centimeters and is provided with a needle to mark & record the daily progress



Samco Medical

Website: www.samcomedical.com,