



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
WEIGHT CUFF SET

Product Code :
SAMCOMEDICAL210003A043



Description :

WEIGHT CUFF SET Strengthens weak muscles through load resisting exercises of both upper and lower limbs. Lead shots filled cuff is double sewn and had long wrap around straps to hold weight cuff securely in Place. A full set consists of 24 color coded cuffs i.e. two each of 150 gm, 300 gm, 500 gm, 750, 1 kg, 1.25 kg, 1.50 kg, 1.75 kg, 2 kg, 2.5 kg, 3 kg and 4 kg



Samco Medical

Website: www.samcomedical.com,