



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
SMOOTH EXERCISE TABLE

Product Code :
SAMCOMEDICAL210003A048



Description :

SMOOTH EXERCISE TABLE Top of 150cm X 85cm laminated board is mounted on rubber tipped tubular legs. Frame finished in oven baked enamel. Height of the top is adjustable from 60cm to 90cm. A cutout enables patient to sit close to the table, providing exercises in horizontal plane for shoulder & elbow with a set of 2 skates



Samco Medical

Website: www.samcomedical.com,