



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
SHOULDER WHEEL

**Product Code :**  
SAMCOMEDICAL210003A052



**Description :**

SHOULDER WHEEL All 16g tube constructed, 100 cm dia wheel revolve smoothly with C.I drum & is fitted with calibrated sensitive resistance mechanism, the resistance is controllable from zero to maximum. A 360 deg. scale is provided on the drum to measure and record the degree of revolution from either direction. Arc of motion can be varied from 26cm to 86cm dia. Fitted with attachment to raise or lower the wheel at the desired height for each patient. Mounted on three laminated hardwood wall boards. Oven baked finish. Packing Details: Size: 43" x 43" x 13". Weight: 50 Kg. each



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),