



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
WHIRLPOOL BATH

**Product Code :**  
SAMCOMEDICAL210005A004



**Description :**

WHIRLPOOL BATH Heavy gauge stainless steel tank sturdily fabricated and bright finished. Provides effective hydro-massage of arms, hips, legs, knees and feet. Mounted on heavy duty 10cm dia four ball bearing castors. Comes complete with removable inside aluminium Seat & Arm rest. Unit is fitted with a motorized turbine ejector and aerator on a spring loaded turbine elevator, digital thermometer, thermostat & immersion heater. Provided with drain pipe for emptying. Inside measurements of tank 90cm long x 50cm wide x 70cm deep. Inside measurements of tank 90cm long x 50cm wide x 70cm deep



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),