



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
HIGH CHAIR

**Product Code :**  
SAMCOMEDICAL210005A006



**Description :**

HIGH CHAIR Used with Whirlpool Baths for Lower Leg & Foot hydro Massage. Broad based tubular steel frame chair is fitted with foam padded, revolving, Seat and Back for height adjustment. Height is adjustable from 80cm to 100cm. Comes with rubber threaded foot step. Fitted with four wheels of 5cm dia. Finished oven baked



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),