



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
ACCESSORIES FOR HYDROTHERAPY TANK

**Product Code :**  
SAMCOMEDICAL210005A007



**Description :**

ACCESSORIES FOR HYDROTHERAPY TANK The optional accessories required for Hydrotherapy tank for conveniently handling of patients. 1. Hoist: An Electric Hoist of 500 kg. capacity for lifting and lowering of patient in the tank. Works on 440 V. three phase A.C. power supply. Hoist is fitted with electric Travelling Trolley that moves the hoist side-ways on the 'I' beam. Comes complete with S.S cross bar with eye hooks for stretcher. 2. Trolley: Painted trolley frame, mounted on four 20cm. dia wheels for carrying stainless steel Stretcher with nylon clothe & belts. NOTE: 'I' beam of 8"x 4" size will provided by the department and got fixed 13 to 15 feet high from floor. Length of beam depends upon the size of the room. Minimum working length 12 feet



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),