



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
HYDROTHERAPY TANK

Product Code :
SAMCOMEDICAL210005A001

Description :

HYDROTHERAPY TANK The compact, rectangular design, full body tank for hydro massage. The hydrotherapy tank with polished interior measures : 200cm long x 100cm wide x 50cm deep from inside. Unit comes with One motorized electric turbine ejector and aerator with turbine carriage and elevator, Digital Thermometer, adjustable Head-Rest fitted with nylon cloth. Tank is fitted with Two water inlets, One overflow and One drain fittings. Tank is finished bright from inside & painted outside



Samco Medical

Website: www.samcomedical.com,