



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Dental Chair Programmable (Conventional Over The Head Delivery)

**Product Code :**  
SAMCOMEDICAL110001A001

**Description :**

Dental Chair Programmable (Conventional Over The Head Delivery)

**Technical Specification :**

Dental Chair Programmable (Conventional Over The Head Delivery) Electrically operated Dental chair unit with Up/Down / Back/ Forward Movement. Feathers touch control on delivery unit & chair backrest. 4 Program mes. 2 ERASABLE GARGLING & AUTO RETURN TO ZERO POSITION. Programmable water unit for spittoon and tumbler through feather touch membrane switches given at front panel. Unit Consisting Followings: Air-Rotor terminal  $\frac{1}{2}$  2 No. Micro - Motor terminal control  $\frac{1}{2}$  1 No. Three way syringe (One at Dr's Assistant Side)-2 No. Two intensity operating light fitted with imported reflector. Arm rest both sides for better patient comfort. Articulated adjustable headrest with facility to convert into pediatric headrest. Stainless Steel instrument Tray. Dr. Stool  $\frac{1}{2}$  Pneumatic. X-Ray viewer  $\frac{1}{2}$  LED based. Seamless Fully laminated imported PU MOULDED upholstery in attractive colors. Low & High Pneumatic suction. Multifunctional foot control. Single piece Ceramic Spittoon top along with Tumbler. 2 Actuators (Noise, Vibration & Jerrck Free). 3 Side Operated Chair



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),