



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Dental Chair Programmable (Conventional Underarm
Delivery Unit)

Product Code :
SAMCOMEDICAL110001A002

Description :

Dental Chair Programmable (Conventional Underarm Delivery Unit)

Technical Specification :

Dental Chair Programmable (Conventional Underarm Delivery Unit) Electrically operated Dental chair unit with Up / Down / Back/ Forward. Auto return to zero position. Feathers touch control on delivery unit & chair backrest. Unit Consisting Followings: Air-Rotor terminal $\frac{1}{2}$ 1 No. Micro - Motor terminal control $\frac{1}{2}$ 1 No. Three way syringe -1 No. Two intensity operating light fitted with imported reflector. Water connection for spittoon and tumbler through electrically operated solenoid valve. Arm rest both sides for better patient comfort. Articulated adjustable headrest with facility to convert into pediatric headrest. Stainless Steel instrument Tray with hanging Straight Tubing. Dr. Stool $\frac{1}{2}$ Pneumatic. X-Ray viewer $\frac{1}{2}$ LED based. Seamless Fully laminated imported upholstery in attractive colors. Multifunctional foot control. Single piece Ceramic Spittoon top. Wide instruments tray. 2 Actuators (Noise, Vibration & Jerrck Free). 3 Side Operated Chair



Samco Medical

Website: www.samcomedical.com,