



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Ultra Sonic Scaler (With 3 Tips)

Product Code :

SAMCOMEDICAL110004A001

Description :

Ultra Sonic Scaler (With 3 Tips)

Technical Specification :

Ultra Sonic Scaler (With 3 Tips) Economic Scaler unit compatible with EMS system. Auto-boost function removes heavy calculus easily. Higher voltage rating components in the circuit reduces its failure rate. Water Cooling System inside keeps working time longer. Top design: all surfaces are easy for clean and maintenance. Autoclavable hand piece can be autoclaved up to 135^o±½C. Frequency: 30KHZ. Power Supply: 110VAC/220VAC, 50/60HZ



Samco Medical

Website: www.samcomedical.com,