



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Dental H.P Oil Free & Noise Less Air Compressor (Fully Imported)

**Product Code :**

SAMCOMEDICAL110002A002

**Description :**

Dental H.P Oil Free & Noise Less Air Compressor (Fully Imported)

**Technical Specification :**

Dental H.P Oil Free & Noise Less Air Compressor (Fully Imported) Heavy-duty air reservoir tank made from thick M.S. sheet with Air & Tank Pressure gauge. Automatic cut off. Switch, safety valve, ON/OFF valve, Drainage valve etc. Moisturizer Filter With Regulator. Tank capacity 50 Lts. V/Hz 200V/50Hz. Easy handling



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),