



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Light Cure Unit With Cordless Digital Lcd Display (LED based)

**Product Code :**

SAMCOMEDICAL110003A002

**Description :**

Light Cure Unit With Cordless Digital Lcd Display (LED based)

**Technical Specification :**

Light Cure Unit With Cordless Digital Lcd Display (LED based) Curing hand piece. Cordless design, LCD display. The special high performance diode mechanism can polymerize an average of 3mm thickness of material within 10 seconds. The smart three setting: Full Power, Soft-up and Impluse Luminosity 1/2 800-1200 mW/cm<sup>2</sup>. LED Lamp: 5 Watt. LED lamp life time up to 100, 000hrs. Wavelength: 450-490nm. Cures in 10 sec./3mm. Additional Lithium-ion Battery



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),