



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ent Unit With Inbuild Recycle Tank

Product Code :
SAMCOMEDICAL110007A001

Description :

Ent Unit With Inbuild Recycle Tank

Technical Specification :

Ent Unit With Inbuild Recycle Tank Standard Equipment LED illuminating light 1 Spray gun (unbent 2 & bend 1) 3 Suction gun 1 Blowing gun 1 Laryngoscope pre-heater 1 Built-in LED light source 1 Luxury thin film viewer 1 Blow-off equipment system (with warning system) 1 Instrument tray 2 Cotton cup 2 Medical bottle 4 Built-in equipment recycle tank 1 Doctor stool 1 Vacuum pump 1 Compressor 1 Option : Imaging system CCD / 3 CCD camera 15 inch LCD monitor Computer (17 inch LCD display / professional endoscope software) Pinter Fiber cable Hernia cold light source Endoscope: ? 2.7 / 4mm, 0ï¿½ / 30ï¿½ otoscope ? 2.7 / 4mm, 0ï¿½ / 30ï¿½ nasoscope ? 8 mm, 70ï¿½ / 90ï¿½ laryngoscope Constant temperature rinsing system LED head light 5 inch intelligent touch controlled



Samco Medical

Website: www.samcomedical.com,