



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Luxury Electric Patient Chair

Product Code :
SAMCOMEDICAL110007A003

Description :

Luxury Electric Patient Chair

Technical Specification :

Luxury Electric Patient Chair The seat can be up & down, backrest can be forward and backward, arms moves together with leg pads, seat can be rotated 180 degree. Voltage : AC220V 50 HZ Power : 230 W (max) and low power consumption, energy saving. Under weight pressure : (minimum) 530mm \times 680mm (maximum) Armrest and backrest connections moving, the backrest can be forwarded 90° into a line Pillow flex : 100mm (folding, detachable) Seat can be rotated 180°



Samco Medical

Website: www.samcomedical.com,