



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
SST (Gel Tubes) (MOQ: 100,000 Pcs Assorted Size)

Product Code :
SAMCOMEDICAL4005A002



Description :

SST (Gel Tubes) (MOQ: 100,000 Pcs Assorted Size) Serum Separating Tubes(SSTs) contain the clot activator and a special gel at the base. The gel forms a stable barrier between the serum and the blood cells. After centrifugation, the serum sample can be used for special diagnostics. Mixing instructions: SSTs should be gently inverted 180° ½ back and forth 5-6 times. Centrifugation instructions: 1300g-2000g for 10 minutes at 25° ½ C room temperature



Samco Medical

Website: www.samcomedical.com,