



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
EDTA K3/K2 Tubes (MOQ: 100,000 Pcs Assorted Size)

Product Code :
SAMCOMEDICAL4005A003



Description :

EDTA K3/K2 Tubes (MOQ: 100,000 Pcs Assorted Size) NET EDTA tubes are used for haematology purposes, where whole blood is required for examination. EDTA functions effectively to cover the calcium ions of the blood, thereby preventing blood coagulation. NET EDTA gel tubes contain EDTA coating and a special gel at the base. The gel forms a stable barrier between the plasma and the blood cells during and after centrifugation. The plasma can be extracted and utilized for examination



Samco Medical

Website: www.samcomedical.com,