



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Sodium Citrate Tubes (MOQ: 100,000 Pcs Assorted Size)

Product Code :
SAMCOMEDICAL4005A004



Description :

Sodium Citrate Tubes (MOQ: 100,000 Pcs Assorted Size) Sodium Citrate Tubes blood collection tubes contain a measured amount of citrate solution and are used for coagulation determinations in citrated plasma. Since the liquid citrate dilutes the blood, it is essential that the tube is completely filled with blood sample so that the dilution is properly carried out. NET Sodim Citrate tubes ensure a constant mixing ratio of one part citrate to nine parts of blood



Samco Medical

Website: www.samcomedical.com,