



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
ESR Vacuum Blood Tubes (MOQ: 100,000 Pcs Assorted Size)

**Product Code :**  
SAMCOMEDICAL4005A006



**Description :**

ESR Vacuum Blood Tubes (MOQ: 100,000 Pcs Assorted Size) ESR is a simple, non-specific screening test utilized for monitoring the presence of inflammation in the body. ESR tubes contain tri-sodium citrate solution and are used for blood sedimentation analysis. Mixing recommendation: ESR tube should be gently inverted 180° ½ back and forth 8-10 times



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),