



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
ESR Vacuum Blood Tubes (MOQ: 100,000 Pcs Assorted Size)

Product Code :
SAMCOMEDICAL4005A006



Description :

ESR Vacuum Blood Tubes (MOQ: 100,000 Pcs Assorted Size) ESR is a simple, non-specific screening test utilized for monitoring the presence of inflammation in the body. ESR tubes contain tri-sodium citrate solution and are used for blood sedimentation analysis. Mixing recommendation: ESR tube should be gently inverted 180 back and forth 8-10 times



Samco Medical

Website: www.samcomedical.com,