



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Multi - Sample Blood Collection Needles With Hub (MOQ: 100,000 Pcs Assorted Size)

**Product Code :**

SAMCOMEDICAL4005A007

**Description :**

Multi - Sample Blood Collection Needles With Hub (MOQ: 100,000 Pcs Assorted Size) The NET Multi -Sample Blood Collection Needles With Hub with optical venipuncture control can make an important contribution to simplifying blood collection routine for phlebotomists. The needle has a transparent plastic hub, which functions as a view window. If venipuncture is successful, the view window is coloured red due to the blood flow



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),