



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Luer Adapter (MOQ: 100,000 Pcs Assorted Size)

Product Code :

SAMCOMEDICAL4005A010

**Description :**

Luer Adapter (MOQ: 100,000 Pcs Assorted Size) The NET Luer adapter is a sterile ,single use device to be used with safety Blood collection sets or luer needles in routine blood collection .It is manufactured from stainless steel and polypropylene ans is fitted with a safety valve at one end.The perforated label not only serves to simplify identification ,but also acts as a seal if integrity



Samco Medical

Website: www.samcomedical.com,