



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Mammography System

Product Code :
SAMCOMEDICAL430001A005



Description :

Mammography System FEATURE : NET Mammography has multi function, very small focus, high resolution, low dosage and is Microprocessor based control. Used for breast screening for iconography examination of mammary diseases and has high sensitivity for mammary benign lumps, nodes, calcified points and early detection of malignant. Finger touch operation panel, easily accessible buttons and foot switch. Intact breast images through streamlined compressor plate. High quality carbon fiber cartridge table with low X-ray absorption. X-ray tube options for the best performance/price ratio: Fixed anode tube Rotational anode tube Italian made lae rotational anode tube American made Varian rotational anode tube The 80KHz high frequency generator designed with advanced IGBT high-frequency inversion technology and quasi-resonant pressure boosting technology. Electrical up/down smooth movement. 360 free C-arm rotation for imaging flexibility



Samco Medical

Website: www.samcomedical.com,